

APPIANFx™ SYSTEM

Biceps Tenodesis Implants

DUAL EXPANDING PEEK ANCHOR

6 mm anchors radially expand to 7.5 mm top, 10 mm bottom for effective tendon-to-bone fixation



Subpectoral

Sliding suture enables gathering of tendon, tension adjustment prior to fixation and additional fixation if desired

Unique Design

- All Inside Arthroscopic Repair
- Minimizes tendon disruption
- Saves O.R. time



Arthroscopic
(Proximal & Suprapectoral)

Cleats facilitate gathering and positioning of tendon

Satisfied Outcomes

- No whip stitching required
- No exteriorizing tendon
- All Inside Repair
- Straightforward delivery and fixation—facilitates a faster procedure

APPIANFx™ SYSTEM


Biceps Tenodesis Implants

DUAL EXPANDING PEEK ANCHOR

Arthroscopic (Proximal & Suprapectoral)



Subpectoral



ORDERING FORMATION

APPIANFx™ Biceps Tenodesis Implant with Delivery Handle	
Model No.	Description
KFX-GTO-160	1 PEEK 6.0mm x 17 mm Biceps Tenodesis Implant w/Delivery Handle

ORDERING FORMATION

APPIANFx™ Biceps Tenodesis Implant with Suture and Delivery Handle	
Model No.	Description
KFX-PSB-150DL	1 PEEK 5.0 mm x 21 mm Biceps Tenodesis Implant w/Force Fiber® Sliding Suture and Delivery Handle
KFX-MST-160	1 PEEK 6.0 mm x 17 mm Biceps Tenodesis Implant w/Force Fiber® Sliding Suture and Delivery Handle



5845 Avenida Encinas, Suite 128
 Carlsbad, CA 92008
 Customer Service: 866.883.8718
 760.444.8844 • Fax: 760.444.8855
www.kfxmedical.com



APPIANFx™ is a trademark of KFX® Medical, Inc

APPIANFx™ is comprised of PEEK-Optima® Polymer from Invibio® Biomaterial Solutions. PEEK-Optima® and Invibio are trademarks of Invibio Ltd. Force Fiber® is a trademark of Teleflex, Inc.
 © 2011 KFX® Medical, All Rights Reserved
 Part No. MK-0067.00 Rev. E, 10-11